# CLINICIAN ARM SUPPORT

## AM0900

## INSTRUCTIONS FOR USE

(Instrucciones de uso / Mode d'emploi / Istruzioni per l'uso / Gebrauchsanweisung / Gebruiksaanwijzing)

Issue Date: 12/2023

Revision: B

EC REP

### EC Representative

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### Manufacturer

Adept Medical Ltd  
2-6 McDonald St,  
Morningside  
Auckland 1025, New Zealand

ADEPT

Medical  
Supporting you

 | [adeptmedicaltraining.com](https://adeptmedicaltraining.com)

[adeptmedical.com](https://adeptmedical.com) | 

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
	<b>A</b>	Arm Support Pad
	<b>B</b>	Forearm Assembly
	<b>C</b>	Tube Assembly
	<b>D</b>	Base Board Assembly

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








	<b>A</b> Arm Support Pad	M3335
	<b>B</b> Forearm Assembly	M9120
	<b>D</b> Base Board Assembly	M9100

### Tube Assembly

### Forearm Assembly (M9120)

## Image Table Compatibility

(Compatibilidad de la mesa de imagen / Compatibilité de la table d'examen / Compatibilità del tavolo per imaging / Kompatibilität Bildgebungstische / Compatibiliteit van tafel)

### Ensure table geometry falls within this specification

- / Asegúrese de que la geometría de la mesa cumpla estas especificaciones
- / Assurez-vous que la géométrie de la table correspond à cette spécification
- / Accertarsi che la geometria del tavolo rientri in questa specifica
- / Stellen Sie sicher, dass die Geometrie des Bildgebungstischs innerhalb dieser Angaben liegt
- / De geometrie van de tafel moet binnen deze specificatie vallen

<b>Table Width:</b> / Anchura de la mesa / Largeur de la table / Larghezza tavolo / Tischbreite / Tafelbreedte	40cm - 70cm
<b>Table Length:</b> / Longitud de la mesa / Longueur de la table / Lunghezza tavolo / Tischlänge / Tafellengte	N/A
<b>Table Angulation:</b> / Angulación de la mesa / Angulation de la table / Inclinazione tavolo / Winkelstellung des Tisches / Tafelhoek	+/- 20° longitudinal & 10° lateral +/- 20° longitudinal y 10° lateral +/- 20° longitudinal et 10° latéral +/- 20° longitudinale e 10° laterale +/- 20° longitudinal und 10° lateral +/- 20° longitudinaal en 10° lateraal
<b>Table Surface:</b> / Superficie de la mesa / Surface de la table / Superficie del tavolo / Tischoberfläche / Tafeloppervlak	Flat Plana Plane Piatto Flach Vlak

# ESSENTIAL INFORMATION

Información básica / Informations essentielles / Informazioni essenziali / Grundlegende Informationen / Essentiële informatie

## 11 Essential Information (EN)

- Information
- Initial Checks & Cleaning
- Disinfection
- Disinfection Instructions
- Approved Cleaners List
- Disposal

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## Information

### Intended Purpose

To support the Clinicians arm during endoscope guided surgery.

### Intended User/Training Requirement

Intended to be used by trained medical professionals.

### Patient Target Group

Patients undergoing endoscopic guided surgery.

### Contraindications

- Not to be used by clinicians weighing over 135 Kg.
- Not to be used with patients weighing under **45 Kg**.

### Warnings and Cautions

- Ensure IFU is read prior to use.
- Ensure Preparatory Cleaning is conducted.
- Ensure Preparatory Product Check is conducted.
- Restless patients could dislodge Clinician Arm Support and need to be monitored.

### Incident Reporting Guidelines

For product complaints and incidents, please complete form on:

[adeptmedical.co.nz/repairs](https://adeptmedical.co.nz/repairs)

## Initial Checks & Cleaning

### Cleaning

Refer to Disinfection Instructions.

### Product Checks

- No visible damage or sharp edges.
- Ensure the **Tube Assembly** is free to move up, down and rotate when the lever is in the unlocked position.
- Ensure the **Tube Assembly** cannot move in any direction when the lever is in the locked position.
- Ensure all rotating joints can be locked / unlocked using thumbscrew.
- **Arm Support Pad**: Ensure no visible damage or tearing.

## Disinfection

### Warnings

- Insufficient cleaning may compromise disinfection process and lead to transmission of pathogens.
- Do not use any disinfectant products or abrasive/corrosive agents which are not on the approved cleaners list.
- Always read manufacturer's instructions and consult the manufacturers MSDS for cleaning and disinfectant products.
- Some disinfectants may cause slight discolouration to the blue material used on some components. This will not affect the strength and the product will remain fit for purpose.

### Limitations on Processing

Discontinue use if:

- any cracks, tears, or breakages are present.
- any sharp edges develop.
- If the locking mechanism no longer holds the **Tube Assembly** in place when locked.

## Disinfection Instructions

### Initial treatment at the point of use

It is important to clean the product once it is removed from the packaging and after each procedure. Ensure all areas including joints, clips, sockets, brackets and levers are thoroughly cleaned after use to remove all contaminant build-up that may be present following a procedure.

### Preparation before cleaning

- Remove the **Arm Support Pad**
- Remove the **Base Board** from beneath the mattress

[adeptmedicaltraining.com/productcleaning/Clinician-Arm-Support](https://adeptmedicaltraining.com/productcleaning/Clinician-Arm-Support)

### Cleaning: Manual

- If any visible residue is present, rinse under water, but avoid submerging the device.
- A soft bristle brush can be used.
- Ensure the device is completely dried prior to disinfection.

### Disinfection

Refer to the Approved Cleaners List for selection of an appropriate disinfecting agent.

- With a cloth dampened by an approved disinfectant, wipe all surfaces.
- With the lever in the unlocked position, move the **Tube Assembly** up and down while wiping to ensure the entire surface has been disinfected
- Ensure the **Forearm Assembly** is fully expanded to ensure all surfaces are disinfected

### Drying

All components should be dried thoroughly before use.

### Maintenance, Inspection and Testing

- Ensure there are no visible damage or sharp edges.
- Ensure the **Tube Assembly** is free to move up and down when the lever is in the unlocked position.
- Ensure the **Tube Assembly** cannot move when the lever is in the locked position.
- **Arm Support Pad**: Ensure no visible damage or tearing.

### Packaging

It is not required to package the device following disinfection.

### Sterilisation

This device should not be subjected to sterilisation processes.

### Storage

Once disinfection is complete and all components are dry the device should be stored in a dry environment.

 *The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.*

## | Approved Cleaners List

### Approved Disinfectant by Brands

#### Liquids:

ORION Laboratories  
70% Isopropyl Alcohol

Jaychem Industries  
2% Chlorhexidine / 70% Alcohol

Betadine  
Povidone Iodine (7.5%)

BODE Chemie GmbH  
Kohrsolin FF

STERIS  
Coverage Spray TB

CaviCide  
Metrex Research

#### Dissolvable:

Du Pont  
Rely+On Virkon Tablets

#### Surface Wipes:

PDI Sani-Cloth Bleach

PDI Sani-Cloth AF3

Reynard Health Supplies  
Surface Disinfectant Wipes

Metrex Research  
CaviWipes

Clinell Wipes - Universal (Green)

Clorox Hydrogen Peroxide Wipes

### Approved Active Disinfectant Ingredients

#### Acid Based:

≤10% Malic Acid CAS 6915-15-7  
≤6% Sulfamic acid CAS 5329-14-6

#### Alcohol Based:

≤5% 2-Butoxyethanol CAS 111-76-2  
≤10% Butyldiglycol CAS 112-34-5  
≤70% Isopropyl Alcohol (Propan-2-ol) CAS 67-63-0  
≤10% Tridecanol CAS 69011-36-5  
<10% Alcohols, C12-14, ethoxylated CAS 68439-50-9  
≤70% Denatured Ethanol CAS 64-17-5

#### Sulphate Based:

≤0.1% PHMB CAS 27083-27-8  
≤55% Potassium Peroxymonosulfate CAS 70693-62-8  
≤3% Potassium Persulfate CAS 7727-21-1

#### Ammonium and Chloride Based:

≤5% Benzalkonium Chloride CAS 68424-85-1  
≤0.28% Benzethonium Chloride CAS 121-54-0  
≤10% Benzyl-C23-18-Alkyl-dimethyl Ammonium Chloride CAS 8001-54-5  
≤2% Chlorhexidine CAS 55-56-1  
≤10% Didecyl Dimethyl Ammonium chloride CAS 7173-51-5  
≤0.5% Quaternary Ammonium Compounds CAS 68956-79-6  
<5% Benzyl-C 12-18 alkyldimethylammonium chlorides CAS 63891-01-5

#### Others:

≤10% Glutaral CAS 111-30-8  
≤0.63% Sodium Hypochlorite CAS 7681-52-9  
≤7.5% Povidone Iodine CAS 25655-41-8  
0.5% Disodium Cocoampho Dipropionate CAS68604-71-7  
≤1.4% Hydrogen Peroxide CAS 7722-84-1

Access the latest Approved Cleaners list via [adeptmedicaltraining.com/downloads](http://adeptmedicaltraining.com/downloads)

## | Disposal

Used product is a biohazard, decontaminate according to instruction provided in the Disinfection Instructions and reach out to the manufacturer on [adeptmedical@adept.co.nz](mailto:adeptmedical@adept.co.nz) for material specification, if required. Dispose according to internal clinic policy taking into consideration local regulations.

[adeptmedical@adept.co.nz](mailto:adeptmedical@adept.co.nz)

### Serious Incidents

Any serious incident which occurs in relation to the device should be reported immediately to:

- The Competent Authority of the applicable Member State.
- The Manufacturer using the Serious Incident Reporting Form found on the Adept Medical Website.

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