PRONE SUPPORT

AM0700

INSTRUCTIONS FOR USE

Issue Date: 03/2026

Revision: H

EC REP

EC Representative

Qualrep Services B.V.
Utrechtseweg 310 Bldg B42
6812 AR Arnhem
Netherlands

Manufacturer

Adept Medical Ltd
2-6 McDonald St,
Morningside
Auckland 1025, New Zealand

 | adeptmedicaltraining.com

ADEPT

Medical
Supporting you

adeptmedical.com | 

Product Registration

Register today to ensure your product is covered by warranty, and for easy access to product services and general support.

www.adeptmedical.com/form/product-registration

Warranty

Service


Support


Contents

Unboxing	4
Assembly	4
Product Use	5
Setup	6
Spare Parts	8
Image Table Verification	9
Essential Information	10


ES Spanish	12	NO Norwegian	34
FR French	14	PL Polish	36
IT Italian	16	PT Portuguese	38
DE German	18	RO Romanian	40
NL Dutch	20	RU Russian	42
HR Croatian	22	SR Serbian	44
CS Czech	24	SK Slovak	46
DA Danish	26	SL Slovenian	48
FI Finnish	28	SV Swedish	50
EL Greek	30	TR Turkish	52
HU Hungarian	32		

Unboxing


	(A)	Face Pad
	(B)	Face Pad Support
	(C)	Adjustable Assembly
	(D)	Base Board


Assembly

4



Disassembly

Product Use

Device Adjustments


Tighten


Loosen


Setup


Adjust the device


Ensure safe clearance to C-arm

Product Removal


Kyphotic Patient - Elevate Torso


Spare Parts

Product Code

	Prone Support	AM0700
--	---------------	--------

Component Codes

	A	Face Pad	M6531
	B	Face Pad Support	M6532
	C	Adjustable Assembly	M6533
	D	Base Board	M6534

Adjustable Assembly (M6533)



Image Table Verification

The product has been risk assessed to operate within the below specifications

CT Bore: ≥ 780 mm
Table Angulation: $\pm 15^\circ$ lateral & longitudinal
Table Surface: Flat / Rectangular

Essential Information

Information

Intended Purpose

To provide a support surface for patient's arm and head in prone position.

Intended User/Training Requirement

Intended to be used by trained medical professionals.

Patient Target Group

Adult male or female patient undergoing procedures in prone position.

Contraindications

- Not to be used with patients over 135 kg.
- Not for paediatric use.
- Not to be used outside of the recommended bore size (refer to page 9).
- Not to be used on patients with extreme Kyphosis or previous spinal injury.
- Not to be used with patients with existing or historical brachial plexus injury.
- Head Assembly not to be used with intubated patient.

Warnings and Cautions

- Ensure IFU is read prior to use.
- Ensure Preparatory Cleaning is conducted.
- Ensure Preparatory Product Check is conducted.
- Ensure patient has adequate mobility for shoulder flexion.
- Depending on patient age, size, weight, mobility, and duration of procedure, the product may not be suitable or may contribute to brachial plexus injury.
- Ensure patients fingers, if over edge, do not impact with Xray camera C-arm.
- Ensure hair is managed in hair net for patients with long hair.
- If comfortable position cannot be achieved or patient is restless, the head assembly can be removed and replaced with a pillow.
- Use of the Product is intended only for contact with intact skin.
- For patients suspected of requiring an emergency transition to the supine position during a procedure, do not use Head Assembly.
- Device is not fixed to table and may dislodge (refer to page 7).
- Before use, ensure patient is on the table and in prone position prior to product placement (refer to page 6).
- Ensure patient's body is appropriately elevated (refer to page 8).
- Check patient for brachial plexus injury regularly.
- Ensure the patient's eyes are clear of the face pad during initial positioning and verify positioning periodically throughout the procedure.
- After use, ensure the product is removed from the table prior to moving patient (refer to page 7).

Incident Reporting Guidelines

For product complaints and incidents, please complete form on www.adeptmedical.co.nz/repairs

Initial Checks & Cleaning

Cleaning

Refer to Disinfection Instructions.

Product Checks

- No break in Carbon Fibre surface.
- No visible damage or sharp edges, e.g. cracks
- Adjustable Assembly: Ensure the assembly is not free to move up and down when lever is in the lock position.
- Adjustable Assembly: Ensure the assembly can move up and down when pushed when lever in the unlocked position.
- Adjustable Assembly: Ensure the assembly is not free to tilt when the screw is locked.
- Adjustable Assembly: Ensure the assembly is free to tilt when the screw is unlocked.
- Ensure Face Pad pin is inserted securely into locating holes.
- Ensure the Adjustable Assembly is fully inserted on the Face Support (refer to page 4).
- Ensure the Adjustable Assembly is fully inserted on the Base Board (refer to page 5).

Disinfection

Warnings

- Insufficient cleaning may compromise disinfection process and lead to transmission of pathogens.
- Do not use any disinfectant products or abrasive/corrosive agents which are not on the approved cleaners list.
- Always read manufacturer's instructions and consult the manufacturers MSDS for cleaning and disinfectant products.
- Do not use any cleaning procedure which exceeds 65° C/149° F.
- Some disinfectants may cause slight discolouration to the soft blue or white material used on some components within the product range. This will not affect the strength and the product will remain fit for purpose.

Limitations on Processing

Discontinue use if:

- Any cracks or breakages are present.
- Any exposed carbon is visible.

Disinfection Instructions

Initial Treatment at the Point of Use

It is important to clean the product once it is removed from the packaging and after each procedure. Ensure all areas including joints, clips, sockets, brackets and levers are thoroughly cleaned after use to remove all contaminant build-up that may be present following a procedure.

Preparation Before Cleaning

- Remove Prone Support from beneath the mattress.
- Remove Adjustable Assembly.
- Remove Face Pad.
- Ensure the Adjustable Assembly is fully extended.

www.adeptmedicaltraining.com/productcleaning/Prone-Support

Cleaning: Manual

- If any visible residue is present, rinse under water but avoid submerging the device.
- A soft bristle brush can be used.
- Ensure the device is completely dried prior to disinfection.

Disinfection

Refer to the Approved Cleaners List for selection of an appropriate disinfecting agent.

Face Pad and Face Support:

- With a cloth dampened by an approved disinfectant, wipe all surface.

Adjustable Assembly:

- With a cloth dampened by an approved disinfectant, wipe all surface, especially the blue lever and screw.

Carbon Fibre Base Board:

- With a cloth dampened by an approved disinfectant, wipe the entire surface, including crevices in which contaminants may build up.

Approved Cleaners List

Approved Disinfectant by Brands

Liquids:

ORION Laboratories - 70% Isopropyl Alcohol
Jaychem Industries - 2% Chlorhexidine / 70% Alcohol
Betadine - Povidone Iodine (7.5%)
BODE Chemie GmbH - Kohrsolin FF
STERIS - Coverage Spray TB
CaviCide - Metrex Research

Dissolvable:

Du Pont - Rely+On Virkon Tablets

Surface Wipes:

PDI Sani-Cloth Bleach
PDI Sani-Cloth AF3
Reynard Health Supplies - Surface Disinfectant Wipes
Metrex Research - CaviWipes
Clinell Wipes - Universal (Green)
Clorox Hydrogen Peroxide Wipes

Approved Active Disinfectant Ingredients

Acid Based:

≤10% Malic Acid CAS 6915-15-7
≤6% Sulfamic acid CAS 5329-14-6

Alcohol Based:

≤5% 2-Butoxyethanol CAS 111-76-2
≤10% Butyldiglycol CAS 112-34-5
≤70% Isopropyl Alcohol (Propan-2-ol) CAS 67-63-0
≤10% Tridecanol CAS 69011-36-5
≤10% Alcohols, C12-14, ethoxylated CAS 68439-50-9
≤70% Denatured Ethanol CAS 64-17-5

Sulphate Based:

≤0.1% PHMB CAS 27083-27-8
≤55% Potassium Peroxymonosulfate CAS 70693-62-8
≤3% Potassium Persulfate CAS 7727-21-1

Ammonium and Chloride Based:

≤5% Benzalkonium Chloride CAS 68424-85-1
≤0.28% Benzethonium Chloride CAS 121-54-0
≤10% Benzyl-C23-18-Alkyl-dimethyl Ammonium Chloride CAS 8001-54-5
≤2% Chlorhexidine CAS 55-56-1
≤10% Didecyl Dimethyl Ammonium chloride CAS 7173-51-5
≤0.5% Quaternary Ammonium Compounds CAS 68956-79-6
≤5% Benzyl-C 12-18 alkyl(dimethylammonium chlorides CAS 63891-01-5

Others:

≤10% Glutaral CAS 111-30-8
≤0.63% Sodium Hypochlorite CAS 7681-52-9
≤7.5% Povidone Iodine CAS 25655-41-8
0.5% Disodium Cocoampho Dipropionate CAS68604-71-7
≤1.4% Hydrogen Peroxide CAS 7722-84-1

Access the latest Approved Cleaners list via www.adeptmedicaltraining.com/downloads

Drying

All components should be dried thoroughly before use.

Maintenance, Inspection and Testing

For all components ensure:

- No exposed Carbon Fibre is visible.
- No visible damage or sharp edges, e.g. cracks.

Face Pad:

- Ensure there is no damage, e.g. tear or pin detachment.

Adjustable Assembly:

- Ensure the assembly is free to move up and down when lever is in unlocked and does not move up and down when lever is locked.
- Ensure the assembly is not free to tilt when the screw is locked and is free to tilt when the screw is unlocked.

Packaging

It is not required to package the device following disinfection.

Sterilisation

This device should not be subjected to sterilisation processes.

Storage

Once disinfection is complete and all components are dry the device should be stored in a dry environment in which the temperature will not exceed 65 °C/149 °F.

Disposal

Used product is a biohazard, decontaminate according to instruction provided in the Disinfection Instructions and reach out to the manufacturer on adeptmedical@adept.co.nz for material specification, if required. Dispose according to internal clinic policy taking into consideration local regulations.

Serious Incidents

Any serious incident which occurs in relation to the device should be reported immediately to:

- The Competent Authority of the applicable Member State.
- The Manufacturer using the Serious Incident Reporting Form found on the Adept Medical Website: www.adeptmedical.com/form/serious-incident-reporting-form

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

Visit the Adept Medical Training Platform for detailed visual instructions related to this device. www.adeptmedicaltraining.com

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