ANTEGRADE IR PLATFORM

AM0440

INSTRUCTIONS FOR USE

Issue Date: 03/2023

Revision: C

EC REP

EC Representative

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Manufacturer

Adept Medical Ltd
2-6 McDonald St,
Morningside
Auckland 1025, New Zealand

ADEPT

Medical
Supporting you

Contents

| | |
|--------------------------|---|
| Unboxing | 4 |
| Assembly | 4 |
| Setup | 4 |
| Product Configurations | 5 |
| Spare Parts | 6 |
| Image Table Verification | 6 |
| Essential Information | 8 |

| | | | |
|----------------|----|-----------------|----|
| ES Spanish | 10 | NO Norwegian | 32 |
| FR French | 12 | PL Polish | 34 |
| IT Italian | 14 | PT Portuguese | 36 |
| DE German | 16 | RO Romanian | 38 |
| NL Dutch | 18 | RU Russian | 40 |
| HR Croatian | 20 | SR Serbian | 42 |
| CS Czech | 22 | SK Slovak | 44 |
| DA Danish | 24 | SL Slovenian | 46 |
| FI Finnish | 26 | SV Swedish | 48 |
| EL Greek | 28 | TR Turkish | 50 |
| HU Hungarian | 30 | | |

Product Registration

Register today to ensure your product is covered by warranty, and for easy access to product services and general support.

www.adeptmedical.com/form/product-registration

Warranty


Service



Support

Unboxing


| | | |
|--|----------|-----------------------------|
| | A | Antegrade Main Platform x 1 |
| | B | Base x 2 |
| | C | Leg x 4 |
| | D | Leg Extension x 2 |



Assembly

Setup



4 Unlock the leg to adjust the platform height

Product Configurations

Drape Support (AM1000)
Not included

Spare Parts

Product Code

| | | |
|---|-----------------------|--------|
|  | Antegrade IR Platform | AM0440 |
|---|-----------------------|--------|

Component Codes


| | | | |
|--|---|-------------------------|-------|
|  | A | Antegrade Main Platform | M3660 |
|  | B | Base | M3662 |
|  | C | Rear Base | M3663 |
|  | D | Leg | M3664 |

Image Table Verification

The product has been risk assessed to operate within the below specifications

Table Angulation: +/- 15° lateral & longitudinal
Table Surface: Flat

Essential Information

Information

Intended Purpose

To provide a work surface over a supine patient for femoral transarterial diagnosis and intervention via Antegrade approach on an Angiography imaging table.

Intended User/Training Requirement

Intended to be used by trained medical professionals.

Patient Target Group

Adult male or female patient requiring diagnostic imaging and image guided therapy.

Contraindications

- Not to be used with patients over 135 kg.
- Not for paediatric use.

Warnings and Cautions

- Ensure IFU is read prior to use.
- Ensure Preparatory Cleaning is conducted.
- Ensure Preparatory Product Check is conducted.
- Ensure all items placed on the Antegrade Main Platform do not exceed a total weight of 10 kg.
- To prevent drape from touching patient face (refer to page 5).

Incident Reporting Guidelines

For product complaints and incidents, please complete form on: www.adeptmedical.co.nz/repairs

Initial Checks & Cleaning

Cleaning

Refer to Disinfection Instructions.

Product Checks

- No exposed Carbon Fibre.
- No visible damage or sharp edges, e.g. cracks.

Leg:

- Rotate all four blue collars to ensure they are free to rotate and provide smooth locking/unlocking.
- Ensure that the Leg is free to move in the unlocked position.
- Ensure that the Leg cannot move up and down when locked.
- Ensure the Legs are securely assembled on the Antegrade Main Platform (refer to page 4).
- Ensure Legs are secured in the Bases (refer to page 4).
- Ensure Legs are locked during procedure (refer to page 5).

Disinfection

Warnings

- Insufficient cleaning may compromise disinfection process and lead to transmission of pathogens.
- Do not use any disinfectant products or abrasive/corrosive agents which are not on the approved cleaners list.
- Always read manufacturer's instructions and consult the manufacturers MSDS for cleaning and disinfectant products.
- Do not use any cleaning procedure which exceeds 65 °C/149 °F.
- Some disinfectants may cause slight discolouration to the blue material used on some components within the product range. This will not affect the strength and the product will remain fit for purpose.

Limitations on Processing

Discontinue use if:

- Any cracks or breakages are present.
- Any exposed carbon fibre is visible.

Disinfection Instructions

Initial treatment at the point of use

It is important to clean the product once it is removed from the packaging and after each procedure. Ensure all areas including joints, clips, sockets, brackets and levers are thoroughly cleaned after use to remove all contaminant build-up that may be present following a procedure.

Preparation before cleaning

- Disassemble Antegrade IR Platform into Leg, Leg Extension, Base and Antegrade Main Platform.
- Remove the product from the table.
- Ensure the Leg is fully extended.

www.adeptmedicaltraining.com/productcleaning/Antegrade-IR-Platform

Cleaning: Manual

- If any visible residue is present, rinse under water, but avoid submerging the device.
- A soft bristle brush can be used.
- Ensure the device is completely dried prior to disinfection.

Disinfection

Refer to the Approved Cleaners List for selection of an appropriate disinfecting agent.

Antegrade Main Platform:

- With a cloth dampened by an approved disinfectant, wipe all surfaces, including the threaded holes on the underside.

Base:

- With a cloth dampened by an approved disinfectant, wipe all surfaces, including leg sockets and blue locking clips.

Leg and Leg Extensions:

- With a cloth dampened by an approved disinfectant, wipe all surfaces.

Approved Cleaners List

Approved Disinfectant by Brands

Liquids:

ORION Laboratories - 70% Isopropyl Alcohol
Jaychem Industries - 2% Chlorhexidine / 70% Alcohol
Betadine - Povidone Iodine (7.5%)
BODE Chemie GmbH - Kohrsolin FF
STERIS - Coverage Spray TB
CaviCide - Metrex Research

Dissolvable:

Du Pont - Rely+On Virkon Tablets

Surface Wipes:

PDI Sani-Cloth Bleach
PDI Sani-Cloth AF3
Reynard Health Supplies - Surface Disinfectant Wipes
Metrex Research - CaviWipes
Clinell Wipes - Universal (Green)
Clorox Hydrogen Peroxide Wipes

Approved Active Disinfectant Ingredients

Acid Based:

≤10% Malic Acid CAS 6915-15-7
≤6% Sulfamic acid CAS 5329-14-6

Alcohol Based:

≤5% 2-Butoxyethanol CAS 111-76-2
≤10% Butyldiglycol CAS 112-34-5
≤70% Isopropyl Alcohol (Propan-2-ol) CAS 67-63-0
≤10% Tridecanol CAS 69011-36-5
<10% Alcohols, C12-14, ethoxylated CAS 68439-50-9
≤70% Denatured Ethanol CAS 64-17-5

Sulphate Based:

≤0.1% PHMB CAS 27083-27-8
≤55% Potassium Peroxymonosulfate CAS 70693-62-8
≤3% Potassium Persulfate CAS 7727-21-1

Ammonium and Chloride Based:

≤5% Benzalkonium Chloride CAS 68424-85-1
≤0.28% Benzethonium Chloride CAS 121-54-0
≤10% Benzyl-C23-18-Alkyl-dimethyl Ammonium Chloride CAS 8001-54-5
≤2% Chlorhexidine CAS 55-56-1
≤10% Didecyl Dimethyl Ammonium chloride CAS 7173-51-5
≤0.5% Quaternary Ammonium Compounds CAS 68956-79-6
<5% Benzyl-C 12-18 alkyl dimethyl ammonium chlorides CAS 63891-01-5

Others:

≤10% Glutaral CAS 111-30-8
≤0.63% Sodium Hypochlorite CAS 7681-52-9
≤7.5% Povidone Iodine CAS 25655-41-8
0.5% Disodium Cocoampho Dipropionate CAS68604-71-7
≤1.4% Hydrogen Peroxide CAS 7722-84-1

Access the latest Approved Cleaners list via www.adeptmedicaltraining.com/downloads

Drying

All components should be dried thoroughly before use.

Maintenance, Inspection and Testing

For all components ensure:

- No visible damage or sharp edges, e.g. cracks.

Antegrade Main Platform:

- No exposed Carbon Fibre is visible.

Base:

- Ensure the locking clips are functional and can retain the Leg.

Leg:

- Ensure thread is intact and can be screwed onto the Antegrade Main Platform.
- Ensure the blue collar rotates smoothly.
- Ensure when blue collar is in the unlocked position the Leg can extend and contract smoothly when force is applied.
- Ensure when blue collar is in the locked position the Leg cannot extend and contract when force is applied.

Packaging

It is not required to package the device following disinfection.

Sterilisation

This device should not be subjected to sterilisation processes.

Storage

Once disinfection is complete and all components are dry the device should be stored in a dry environment.

Disposal

Used product is a biohazard, decontaminate according to instruction provided in the Disinfection Instructions and reach out to the manufacturer on adeptmedical@adept.co.nz for material specification, if required. Dispose according to internal clinic policy taking into consideration local regulations.

Serious Incidents

Any serious incident which occurs in relation to the device should be reported immediately to:

- The Competent Authority of the applicable Member State.
- The Manufacturer using the Serious Incident Reporting Form found on the Adept Medical Website: www.adeptmedical.com/form/serious-incident-reporting-form

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

Visit the Adept Medical Training Platform for detailed visual instructions related to this device. www.adeptmedicaltraining.com

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