
MR Safe

OVERHEAD ARM SUPPORT

AM4000

INSTRUCTIONS FOR USE

Issue Date: 08/2024

Revision: E

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Manufacturer
Adept Medical Ltd
2-6 McDonald St,
Morningside
Auckland 1025,
New Zealand

Product Registration

Register today to ensure your product is covered by warranty, and for easy access to product services and general support.

www.adeptmedical.com/form/product-registration

Warranty

Service

Support

Product Code Warning


- Overhead Arm Support with an **AM2000** Product Code contains an inner metal spring and IS NOT MRI compatible.
- Overhead Arm Support with an **AM4000** product code has passed MRI compatibility testing and is safe to use as a patient support in an MR imaging centre with a bore size 70cm and upward.

Contents


Unboxing	4
Assembly	4
Setup	4
Spare Parts	6
Image Table Verification	7
Essential Information	8

ES Spanish	10	NO Norwegian	32
FR French	12	PL Polish	34
IT Italian	14	PT Portuguese	36
DE German	16	RO Romanian	38
NL Dutch	18	RU Russian	40
HR Croatian	20	SR Serbian	42
CS Czech	22	SK Slovak	44
DA Danish	24	SL Slovenian	46
FI Finnish	26	SV Swedish	48
EL Greek	28	TR Turkish	50
HU Hungarian	30		



Unboxing


Assembly


Setup


Height Adjustment


Strapping (Optional)

Strap Placement Options

Padding (Optional)

CT/MRI Safety Clearance



Image Table Verification

The product has been risk assessed to operate within the below specifications

Table Angulation: +/- 15° lateral & longitudinal
Table Surface: Flat / Curved
CT/MRI Bore: ≥70cm

Spare Parts

Product Code

	Overhead Arm Support	AM4000
--	----------------------	--------

Component Codes

	A Armrest Wing MR Safe	M5504
	B Arm Strap x 10	AM0620
	C Locking Leg Assembly MR Safe	M5505
	D Daggerboard MR Safe	M5506

Essential Information

Information

Intended Purpose

To support the patient's arms overhead with controlled shoulder flexion for imaging. Intended for bores 70cm and upward.

Intended User/Training Requirement

Intended to be used by trained medical professionals.

Patient Target Group

Adult male or female patient requiring diagnostic imaging and image guided therapy.

Contraindications

- Not to be used with patients over 135 kg.
- Not for paediatric use.
- Not intended for bore size less than 70cm.
- Not to be used with patients with existing or historical brachial plexus injury.

Warnings and Cautions

- Ensure IFU is read prior to use.
- Ensure Preparatory Cleaning is conducted.
- Ensure Preparatory Product Check is conducted.
- Ensure MRI safety checks are conducted.
- Ensure patient and device will fit clear of the bore.
- Do not use with patients over 135 kg.
- Not for paediatric use.
- Ensure general patient insulation protocol is followed for MRI.
- Padding recommended on support surfaces (refer to page 6).
- For CT and MRI, ensure there is adequate patient clearance prior to use (refer to page 6).
- Depending on patient age, weight, mobility, and duration of procedure, the product may contribute to brachial plexus injury.
- Check patient for brachial plexus injury regularly.

Incident Reporting Guidelines

For product complaints and incidents, please complete form on: www.adeptmedical.co.nz/repairs

Initial Checks & Cleaning

Cleaning

Refer to Disinfection Instructions.


Product Checks

- No visible damage or sharp edges, e.g. cracks.
- Arm Strap: not damaged, e.g. tearing.
- Locking Leg Assembly: ensure that the Locking Leg Assembly is not free to move in the locked position.
- Locking Leg Assembly: Ensure it can move up and down when pushed in the unlocked position.
- Ensure the Locking Leg Assembly is fully inserted into the Daggerboard and the release clip is fully engaged (refer to page 4).
- Ensure the Armrest Wing is fully inserted onto the Locking Leg Assembly (refer to page 4).


MRI Safety Checks

Ensure each component below has the 'MR Safe' product sticker and engraving in the below locations:


- Armrest Wing:


- Locking Leg Assembly:

- Daggerboard:

- Safety Clearance: ≥70cm

Patient Management

Ensure insulation protocol is followed according to individual facility guidelines. This includes ensuring sufficient insulation is placed to eliminate skin to skin conductive loops and that the patient's skin does not contact the bore of the magnet, when using the device.

Disinfection

Warnings

- Insufficient cleaning may compromise disinfection process and lead to transmission of pathogens.
- Do not use any disinfectant products or abrasive/corrosive agents which are not on the approved cleaners list.
- Always read manufacturer's instructions and consult the manufacturers MSDS for cleaning and disinfectant products.
- Some disinfectants may cause slight discolouration to the soft blue material used on some components within the product range. This will not affect the strength and the product will remain fit for purpose.

Limitations on Processing

Discontinue use if:

- Any cracks or breakages are present.

Disinfection Instructions

Initial treatment at the point of use

It is important to clean the product once it is removed from the packaging and after each procedure. Ensure all areas including joints, clips, sockets, brackets and levers are thoroughly cleaned after use to remove all contaminant build-up that may be present following a procedure.

Preparation before cleaning

- Remove Overhead Arm Support from beneath the mattress.
- Remove Arm Straps from Armrest Wing.
- Disassemble Overhead Arm Support into Armrest Wing, Locking Leg Assembly and Daggerboard.
- Ensure the Locking Leg Assembly is fully extended.

Cleaning: Manual

- If any visible residue is present, rinse under water, but avoid submerging the device.
- A soft bristle brush can be used.
- Ensure the device is completely dried prior to disinfection.

Disinfection

Refer to the Approved Cleaners List for selection of an appropriate disinfecting agent.

Armrest Wing

- With a cloth dampened by an approved disinfectant, wipe all surfaces, including inside the Arm Strap slots.

Daggerboard

- With a cloth dampened by an approved disinfectant, wipe all surfaces, including inside of the Locking Leg Assembly slot.

Locking Leg Assembly

- With a cloth dampened by an approved disinfectant, wipe all surfaces, especially the blue lever.

Arm Strap

- With a cloth dampened by an approved disinfectant, wipe all surfaces.

www.adeptmedicaltraining.com/productcleaning/Overhead-Arm-Support

Approved Cleaners List

Approved Disinfectant by Brands

Liquids:

ORION Laboratories - 70% Isopropyl Alcohol
Jaychem Industries - 2% Chlorhexidine / 70% Alcohol
Betadine - Povidone Iodine (7.5%)
BODE Chemie GmbH - Kohrsolin FF
STERIS - Coverage Spray TB
CaviCide - Metrex Research

Dissolvable:

Du Pont - Rely+On Virkon Tablets

Surface Wipes:

PDI Sani-Cloth Bleach
PDI Sani-Cloth AF3
Reynard Health Supplies - Surface Disinfectant Wipes
Metrex Research - CaviWipes
Clinell Wipes - Universal (Green)
Clorox Hydrogen Peroxide Wipes

Approved Active Disinfectant Ingredients

Acid Based:

≤10% Malic Acid CAS 6915-15-7
≤6% Sulfamic acid CAS 5329-14-6

Alcohol Based:

≤5% 2-Butoxyethanol CAS 111-76-2
≤10% Butyldiglycol CAS 112-34-5
≤70% Isopropyl Alcohol (Propan-2-ol) CAS 67-63-0
≤10% Tridecanol CAS 69011-36-5
≤10% Alcohols, C12-14, ethoxylated CAS 68439-50-9
≤70% Denatured Ethanol CAS 64-17-5

Sulphate Based:

≤0.1% PHMB CAS 27083-27-8
≤55% Potassium Peroxymonosulfate CAS 70693-62-8
≤3% Potassium Persulfate CAS 7727-21-1

Ammonium and Chloride Based:

≤5% Benzalkonium Chloride CAS 68424-85-1
≤0.28% Benzethonium Chloride CAS 121-54-0
≤10% Benzyl-C23-18-Alkyl-dimethyl Ammonium Chloride CAS 8001-54-5
≤2% Chlorhexidine CAS 55-56-1
≤10% Didecyl Dimethyl Ammonium chloride CAS 7173-51-5
≤0.5% Quaternary Ammonium Compounds CAS 68956-79-6
≤5% Benzyl-C 12-18 alkylidimethylammonium chlorides CAS 63891-01-5

Others:

≤10% Glutaral CAS 111-30-8
≤0.63% Sodium Hypochlorite CAS 7681-52-9
≤7.5% Povidone Iodine CAS 25655-41-8
0.5% Disodium Cocoampho Dipropionate CAS68604-71-7
≤1.4% Hydrogen Peroxide CAS 7722-84-1

Access the latest Approved Cleaners list via www.adeptmedicaltraining.com/downloads

Drying

All components should be dried thoroughly before use.

Maintenance, Inspection and Testing

For all components ensure:

- No visible damage or sharp edges, e.g. cracks.

Locking Leg Assembly:

- Ensure blue lever functions smoothly.
- Ensure when lever is in the unlocked position the Locking Leg Assembly can extend and contract smoothly when force is applied.
- Ensure when lever is in the locked position the Locking Leg Assembly cannot contract or extend when force is applied.
- Ensure smooth engagement with daggerboard during assembly.
- Ensure tight fit with Armrest Wing during assembly.

Packaging

It is not required to package the device following disinfection.

Sterilisation

This device should not be subjected to sterilisation processes.

Storage

Once disinfection is complete and all components are dry the device should be stored in a dry environment.

Disposal

Used product is a biohazard, decontaminate according to instruction provided in the Disinfection Instructions and reach out to the manufacturer on adeptmedical@adept.co.nz for material specification, if required. Dispose according to internal clinic policy taking into consideration local regulations.

Serious Incidents

Any serious incident which occurs in relation to the device should be reported immediately to:

- The Competent Authority of the applicable Member State.
- The Manufacturer using the Serious Incident Reporting Form found on the Adept Medical Website: www.adeptmedical.com/form/serious-incident-reporting-form

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

Visit the Adept Medical Training Platform for detailed visual instructions related to this device. www.adeptmedicaltraining.com

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